



 **smith&nephew**
EXU-DRY[◊]
Anti Shear Highly Absorbent

Anti-shear. Providing your patients
extra layers of comfort and protection.

Anti-shear. Providing your patients extra layers of comfort and protection.

- Non-adherent
- Semi-Permeable, non-occlusive
- Anti-Shear, highly absorbent
- Soft and pliable
- One piece, multi-layer

EXU-DRY is a one-piece, multi-layer, highly absorbent, non-adherent wound dressing. Its construction incorporates a non-adherent wound contact layer which makes dressing changes less painful and traumatic; an anti-shear layer which reduces friction and shearing to the wound; absorbent inner layers which wick exudate away from the wound bed; and a permeable, non-occlusive outer layer which allows the wound to breathe.

The sealed, multi-layer construction of EXU-DRY means nurses do not need to apply their own series of layers, therefore reducing nursing time. In addition, its high absorbency reduces the number of dressing changes required.

Advantages

Highly Absorbent

EXU-DRY is highly absorbent and holds many times its weight in fluid. This reduces the number of dressing changes required and therefore reduces nursing time. The inner absorbent layers wick exudate away from the wound, reducing the risk of maceration of the surrounding skin.^{1,2}

Non-adherent

The wound contact layer prevents adherence to the wound and reduces trauma to the wound bed when the dressing is removed.²

Anti-shear

EXU-DRY has an anti-shear layer which reduces friction and shearing during patient movement. This allows the wound bed to heal undisturbed.²

One-piece, Multi-layer

EXU-DRY incorporates many layers into one dressing, and simply requires securing in place, which enables easier and quicker dressing changes for the nurse and patient. The sealed construction of EXU-DRY eliminates frayed edges and will not leave debris in the wound.

Semi-permeable, Non-occlusive Outer Layer

The outer layer permits the wound to breathe, thus preventing maceration and breakdown of the surrounding intact skin.²

Soft and Pliable

EXU-DRY is soft and pliable to ensure effective contact with the wound and provide excellent patient comfort.

Ordering information

Code	Product description	Items per unit
5999003	7.5cm disc; medium	Box/50
5999034	7.5cm x 10cm; full	Box/50*
5999M04	10cm x 15cm; medium	Box/10***
5999004120	10cm x 15cm; full	Box/10
5999006	15cm x 23cm; full	Box/12
5999M06	15cm x 23cm; medium	1 Dressing***
5999009	23cm x 38cm; full	Box/10
5999018	38cm x 46cm; full	Box/30**
5999M18	38cm x 46cm; medium	1 Dressing***
5999024	38cm x 61cm; full	Box/30*
5999M24	38cm x 61cm; medium	1 Dressing***
5999028	50cm x 71cm; full	Box/20*
5999M28	50cm x 71cm; medium	Box/30*
Slit Tube Dressings		
5999003S	8cm disc; medium	Box/50*
5999PTM	5cm x 8cm; full	Box/50*
5999034S	8cm x 10cm; full	Box/50*
5999101	10cm x 15cm; full	1 Dressing***
Incision Dressings		
5999INC9	8cm x 22cm; full	Box/50*
Specialty Dressings		
5999LPG	35cm x 35cm; hand dressing (glove shaped); heavy	Box/20
5999LPL	87cm x 95cm; leg dressing (shaped); heavy	Box/20*
5999M36	61cm x 91cm; non-permeable; heavy	Box/15*
5999L72	91cm x 183cm; non-permeable	1 Dressing***
5999L73	91cm x 183cm; non-permeable with straps	Box/15*

*Not available in New Zealand

**Available for purchase as a single unit in New Zealand

***Not available in Australia

References

1. Edwards J. Use of Exu-Dry in the management of a variety of exuding wounds. *Brit J Nurs* 2001; 10(12):815-818.
2. Brown-Eltris M. Considering dressing changes. *Ostomy Wound Man* 1994; 40:46-51.

Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999 F +61 2 9857 3900
www.smith-nephew.com/australia
Customer Service
T 13 13 60 F 1800 671 000

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840 F +64 9 820 2841
www.smith-nephew.com/new-zealand
Customer Service
T 0800 657 799 (Surgical) T 0800 807 663 (Wound)



®Trademark of Smith & Nephew
SNI2867 (11/16)